

BETTER WATER, INC[®]

BICARB SYSTEMS - MEDICAL-DEVICE

BICARB SYSTEMS



Central delivery bicarbs are used in the mixing and delivery of Sodium Bicarbonate Solution. Sodium Bicarbonate is part of the final dialysate solution that is given to the individual dialysis patients.

The central delivery bicarb system is designed with patient safety, ease of installation and versatility in mind.

The unit is pre-plumbed. The only necessary plumbing required for installation is a dialysis water hookup and a bicarb loop feed and return. All that is necessary to connect the electrical part of the system, is to simply plug in the accompanying 115 vac plug into the specified receptacle.

The bicarb system carries a 510k certificate from the FDA. Each unit is individually tested prior to shipping. Equipment includes enclosed UL 508A , industrial panels.

Standard Features

- High level shut off
- Low tank alarm
- Mix and/or feed from either tank (Dual Systems)
- Polyethylene tank and stand



Solutions. Service. Support
www.betterwater.com

Equipment changes are at the sole discretion of BWI. Brochure photographs & specs may not reflect recent changes.

BICARB SYSTEM

MODELS AND SPECIFICATIONS

Model	EQB-61H-100	EQB-62H-100	EQB-11H-100	EQB-12H-100
Total Capacity	60 gal-	120 gal-	100 gal-	200 gal-
# of	1	2	1	2
Tank and Stand Material	Polyethylene			
Control Voltage	24			
Overall Dimensions	2'3" x 3'6"	2'3" x 5'8"	3'3" x 4'0"	3'3" x 6'6"
Est. Operating Weight	810 lbs	1349	1240	2065
Electrical Require-	110 vac/20 amps/60 Hz			
Bicarb Loop Connections	¾" pvc sch 80 Female Socket, solvent weld			
Pure Water Loop Connections	1" pvc sch 80 Female Socket, solvent			
Additional Bicarb Options Available				

BWI carries over 2000 parts, accessories and consumables in stock, and ready to ship including:

Carbon Filters (provides 10 Minute Empty Bed Contact Time)
 Duplex & Simplex Softeners
 Storage Tanks
 Softener
 Brass pressure regulator
 Blend valve
 City boost pump
 Thermometer
 Floor valve boxes
 Remote alarms
 & More!



Central Water Systems Central Delivery Systems Reverse Osmosis Units Portable RO Units

Pre-Treatment Racks Post Treatment Racks D.I. Exchange Accessories Consumables

Hemodialysis Medical Pharmaceutical Laboratories Commercial Industrial

Better Water, Inc
 698 Swan Drive Smyrna, TN 37167
 Phone: 615.355.6063 Fax: 615.355.6065
www.betterwater.com