

BETTER WATER, INC[®]

PRE-TREATMENT RACK - MEDICAL DEVICE



EcoSmart PRE-TREATMENT

Smart for you, smart for the planet!

MOBILE PRE-TREATMENT RACK

Better Water, Inc's uniquely designed Pre-Treatment Rack, provides another level of mobility to your water treatment system. Instead of pre-treatment components being wall mounted, they mount on a single frame, mobile unit.

Instead of being assembled on site, the pre-treatment components are manufactured in a controlled environment. This allows for unparalleled quality and consistency for the pre-treatment portion of the water system. Pre-plumbed and pre-wired, the unit decreases installation time. Equipment includes enclosed, UL 508A industrial panels. The new energy efficient design **provides up to a 50% energy savings.**

Standard Features

- Blend Valve
- Blend Valve By-Pass
- Multi-Cartridge Housing
- City Boost Pump
- 2 Auxiliary Ports
- Copper Piping
- Gauges
- Control Box
- Rugged metal frame construction
- Baked on powder coat finish
- Standard alarms
- Mobile - locking wheels



Energy Savings



Solutions. Service. Support
www.betterwater.com

Equipment changes are at the sole discretion of BWI. Brochure photographs & specs may not reflect recent changes.

PRE-TREAT RACK

MODELS AND SPECIFICATIONS					
	EQPRE-410-CE-	EQPRE-410-CK-	EQPRE-515-CK-	EQPRE-715-CK-	EQPRE-715-CL-
CAPACITY	18gpm	18gpm	22gpm	25gpm	25gpm
PIPING SIZE	1"	1"	1 1/2"	1 1/2"	1 1/2"
BLEND VALVE	SYMMONS 5-400	SYMMONS 5-400	SYMMONS 5-500	SYMMONS 5-700	SYMMONS 5-700
BLEND VALVE	YES	YES	YES	YES	YES
CITY BOOSTER	2hp, 1ph	2hp, 3ph	2hp, 3ph	2hp, 3ph	3hp, 3ph
GAUGES	2 1/2", all Stainless Steel, 1.5% Accuracy, Dry, 0-100psi. Incoming Hot, Incoming Cold, Post Blend Valve, Post Multi-Cartridge Housing, Post City Booster Pump.				
FILTRATION	All Stainless Steel, Multi-Cartridge Housing, Containing 4) 20", 5µ Filters. (other µ rating filters available)				
CONTROL BOX	Control Box on face of unit, utilizing relay logic controls				
SPACE REQUIRE-	24"W x 60"D	24"W x 60"D	24"W x 60"D	24"W x 60"D	24"W x 60"D
POWER	230vac, 1ph, 60Hz,	208vac, 3ph, 60Hz,	208vac, 3ph, 60Hz,	208vac, 3ph, 60Hz,	208vac, 3ph, 60Hz,
INCOMING WATER REQUIREMENTS	1" Hot, 1" Cold, Dedicated Copper lines	1" Hot, 1" Cold, Dedicated Copper lines	1 1/2" Hot, 1 1/2" Cold, Dedicated Copper lines	1 1/2" Hot, 1 1/2" Cold, Dedicated Copper lines	1 1/2" Hot, 1 1/2" Cold, Dedicated Copper lines
DRAIN REQUIREMENTS	30 gpm (Floor Drain) or consult factory for sizing	30 gpm (Floor Drain) or consult factory for sizing	30 gpm (Floor Drain) or consult factory for sizing	30 gpm (Floor Drain) or consult factory for sizing	30 gpm (Floor Drain) or consult factory for sizing

**Optional, Dual City Booster Pumps, Sediment Filtration , and provisions for Break Tank Applications available for all models.

Additional models are available. Please contact BWI for information.

BWI carries over 2000 parts, accessories and consumables in stock, and ready to ship including:

Carbon Filters (provides 10 Minute Empty Bed Contact Time)

Duplex and Simplex Softeners

Storage Tanks

Brass pressure regulator

Blend valve

City boost pump

Thermometer

Floor valve boxes

Remote alarms

& More!



Central Water Systems Central Delivery Systems Reverse Osmosis Units Portable RO Units

Pre-Treatment Racks Post Treatment Racks D.I. Exchange Accessories Consumables

Hemodialysis Medical Pharmaceutical Laboratories Commercial Industrial

Better Water, Inc
 698 Swan Drive Smyrna, TN 37167
 Phone: 615.355.6063 Fax: 615.355.6065
 www.betterwater.com