

BETTER WATER LLC

Floor Valve Boxes

Medical Device

FLOOR BOX



For effective containment of dialysis machine connections, Better Water LLC has designed the Floor Valve Box. This component serves as an easily accessible container which conveniently houses all your dialysis machine connections. With the BWI Floor Valve Box your dialysis work area is kept free from leakage and corrosion.

The Better Water LLC Floor Valve Box is made of light-weight thermoplastic. Thermoplastic is corrosion-free, durable and easily cleaned. The Better Water Floor Valve Box is custom built for your application and can accommodate up to four loops for primary and secondary concentrate, bicarbonate and pure water. For easy drain connections, the unit is equipped with a drain tail piece.

Standard Features

- Easy access to machine connections
- Made of light-weight Thermoplastic
- Effectively contains leakage
- Corrosion Free
- Connections are quick connect / disconnect
- Built in check valves to minimize spillage
- Concentrate and Bicarbonate regulators are included
- Tubing, fittings and valves are included for connection to loops


Solutions. Service. Support.

Better Water LLC 698 Swan Dr Smyrna, Tn 37167 (615) 355-6063 and at www.betterwater.com

Equipment changes are at the sole discretion of BW,LLC. Brochure photographs & specs may not reflect recent changes



Floor Valve Box



BW LLC carries over 2000 parts, accessories and consumables in stock, and ready to ship including:

- Carbon Filters (provides 10 Minute Empty Bed Contact Time)
- Duplex & Simplex Softener
- Storage Tanks
- Brass pressure regulator
- Blend valve
- City boost pump
- Thermometer
- Floor valve boxes
- Remote alarms
- & More!



Central Water Systems	Central Delivery Systems	Reverse Osmosis Units	Portable RO Units
Pre-Treatment Racks	Post Treatment Racks	D.I. Exchange Accessories	Consumables
Hemodialysis	Medical	Pharmaceutical	Laboratories
		Commercial	Industrial