

BETTER WATER LLC


Portable R.O. Carts

Portable R.O. Carts

Better Water LLC has your solution for making your portable RO system truly portable. Our Short and Extended carts allow for easy moving of the RO and dialysis machine and the supporting carbon and softener tanks from location to location. Strong, light weight and silent, these carts are the perfect long lasting accessory allowing easy one person movement and when it's not being moved, you'll rest assured it will be where you left it due to its locking casters.

With
Stainless
Steel
Frame!

Short cart



Extended cart

Standard Features

- Rugged, long lasting Stainless Steel frame
- Bead blasted
- Rust inhibited
- Silent non-marking wheels
- Two locking casters
- Short carts for the portable RO and pretreatment
- Extended carts will support the RO, pretreatment AND the Dialysis machine.



Solutions. Service. Support.

Better Water LLC 698 Swan Dr Smyrna, Tn 37167 (615) 355-6063 and at www.betterwater.com

Equipment changes are at the sole discretion of BW,LLC. Brochure photographs & specs may not reflect recent changes

Rev. June 2011