

BETTER WATER LLC

Post Treatment Assemblies

Post Treatment Assembly

Sized accordingly

Better Water LLC has the solution to all your post Reverse osmosis and final filtration needs. The Post Treatment assemblies accommodate re-pressurization Pumps, D.I. connections, D.I. resin trap, final filtration, and all the necessary controls and safety features. Each assembly has the same features and safety considerations of Better Water's standard water system, and can be pre-fabricated on a single frame or assembled on site as a wall mount. This allows for unparalleled flexibility. One of the most important parts of the water treatment system-Final Filtration

The assemblies come assembled or ready for installation. Simply plumb to a storage tank, the R.O. water loop feed and to the drain. This makes for a perfect fit in any new water treatment system installation or retrofitting to an existing facility.



Shown with 4x40 Ultra Filters

Standard Features

- Dual alternating Repress pumps
- Flow Switch
- Ultra filter assembly
- .03 or 4x40
- D.I. provisions
- Divert to drain
- Final water quality monitor (conductivity)



Shown with .03 Filters



Solutions. Service. Support.

Better Water LLC 698 Swan Dr Smyrna, Tn 37167 (615) 355-6063 and at www.betterwater.com

Equipment changes are at the sole discretion of BW,LLC. Brochure photographs & specs may not reflect recent changes

Post Treatment Assemblies



Models and Specifications		
Model (examples only)	EQPST-255-RM-GJ	EQPST-303-RM-GJ
Capacity GPM	12	18
# of Filters/Membranes	3	3
Filters/Membranes	4x40 spiral wound	20" .03
Operating Weight	425 lbs	290 lbs
Overall Dimensions	25" w x 45" d x 72" h	
Power Requirements	208vac, 3ph, 60Hz	
Number of Pumps (Auto-Alternating)	2	2
Incoming water requirements	Reverse Osmosis Water	
Gauges	All stainless steel construction, Pre and Post gauges.	
Other standard features	Sample ports, resin trap, flow switch ,compatible with 1% Peracetic acid and Bleach for disinfecting	

BW LLC carries over 2000 parts, accessories and consumables in stock, and ready to ship including:

- Carbon Filters (provides 10 Minute Empty Bed Contact Time)
- Duplex & Simplex Softener
- Storage Tanks
- Brass pressure regulator
- Blend valve
- City boost pump
- Thermometer
- Floor valve boxes
- Remote alarms
- & More!



Central Water Systems	Central Delivery Systems	Reverse Osmosis Units	Portable RO Units
Pre-Treatment Racks	Post Treatment Racks	D.I. Exchange Accessories	Consumables
Hemodialysis	Medical	Pharmaceutical	Laboratories
		Commercial	Industrial