

BETTER WATER LLC

Pre-Treatment An Overview



Pre-Treatment Rack



Softener

PRE TREATMENT is one of the basic components in a water system. The blending valve through the softener is part of the pre - treatment section. The entire water system's operation depends on proper pretreatment.

Proper pre - treatment can lengthen the life of the R.O. membranes, improve the quality of water produced, and reduce the amount of maintenance and cleaning a system may require.

Using our Decision Making Process (DMP), BW LLC reviews the feed water analysis along with the clinical demands for RO water and determines the size and types of pre - treatment required to give your water system optimal performance.



Pre-Treatment Wall Mount Header

Examples of Pre-Treatment Components

- Blend Valve
- City Booster Pumps
- Carbons
- Softeners
- Sediment Filtration

Accurate pre - treatment, scheduled feed water testing, and product water monitoring, along with proper cleaning and disinfecting, are preventive measures that help to ensure the quality of the product water and extend the life of your RO system.

Solutions. Service. Support.

Better Water LLC 698 Swan Dr Smyrna, Tn 37167 (615) 355-6063 and at www.betterwater.com

Equipment changes are at the sole discretion of BW,LLC. Brochure photographs & specs may not reflect recent changes

Rev. June 2011

