

# BETTER WATER, INC<sup>®</sup>

MEDIPAC RO SYSTEM - MEDICAL DEVICE

## MediPac RO System



The MediPAC water treatment system is pre-engineered and pre-assembled water system designed with easy installation and versatility in mind. The MediPAC comes on casters and is entirely self-contained. The unit is pre-plumbed. The only necessary plumbing needed for installation is a hot and cold water supply, loop feed and return, plumbed reservoir and to connect the system drain. The wiring is also completed at the factory, all that is required to hook-up the electrical part of the system is to simply plug the accompanying 115 vac power cord into the specified receptacles. A remote alarm box is included with the system.

The MediPac System carries a 510k certificate from the FDA and is designed to meet or exceed AAMI RD: 62 on RO quality water for dialysis. Each unit is individually tested prior to shipping.



## Standard Features

- Tap water boost pump and control (if needed)
- Blend valve complete with temperature gauge and cold water by-pass
- 0-100 PSI pressure gauges pre and post every piece of equipment
- Automatic back washing carbon filters
- Automatic regenerating water softener (if needed)
- Pre-R.O. 5 micron filtration
- Interlock circuit with indicator lights for the automatic back washing equipment
- Tank feed systems include FDA grade polyethylene or polypropylene R.O. water storage tank, complete with air filter, low level alarm and controls
- Repressurization pump with control box and flow switch (if needed)



**Solutions. Service. Support**  
**[www.betterwater.com](http://www.betterwater.com)**

Equipment changes are at the sole discretion of BWI. Brochure photographs & specs may not reflect recent changes.

# MEDIPAC SYSTEM

## MEDIPAC WATER TREATMENT SYSTEMS FOR HEMODIALYSIS

Model	Tank or Direct Feed	Blend Valve	Permeate Flow @77°F (+/-15%)	Softener Grains Capacity	Softener Size	Total Cu. Ft. of Carbon	Size of Carbon Filters	Capacity of Storage Tank
HD-2200TF	TANK	YES	1.5 GPM	60,000	12" x 52"	4.0	12" x 52"	60 GAL
HD-2200DF	DIRECT	YES	1.5 GPM	60,000	12" x 52"	4.0	12" x 52"	NA
HD-2880TF	TANK	YES	2.0 GPM	60,000	12" x 52"	6.0	14" x 65"	60 GAL
HD-2880DF	DIRECT	YES	2.0 GPM	60,000	12" x 52"	6.0	14" x 65"	NA

*Custom configurations available*

***BWI carries over 2000 parts, accessories and consumables in stock, and ready to ship including***

Carbon Filters (provides 10 Minute Empty Bed Contact Time)  
 Storage Tanks  
 Duplex & Simplex Softener  
 Brass pressure regulator  
 Blend valve  
 City boost pump  
 Thermometer  
 Floor valve boxes  
 Remote alarms  
 & More!



Central Water Systems    Central Delivery Systems    Reverse Osmosis Units    Portable RO Units

Pre-Treatment Racks    Post Treatment Racks    D.I. Exchange    Accessories    Consumables

Hemodialysis    Medical    Pharmaceutical    Laboratories    Commercial    Industrial

Better Water, Inc  
 698 Swan Drive Smyrna, TN 37167  
 Phone: 615.355.6063    Fax: 615.355.6065  
[www.betterwater.com](http://www.betterwater.com)