

BETTER WATER, INC[®]

PRE-TREATMENT - An Overview

Pre-Treatment Rack



PRE TREATMENT is one of the basic components in a water system. The blending valve through the softener is part of the pre-treatment section. The entire water system's operation depends on proper pretreatment.

Proper pre-treatment can lengthen the life of the R.O. membranes, improve the quality of water produced, and reduce the amount of maintenance and cleaning a system may require.

Using our Decision Making Process (DMP's), BWI reviews the feed water analysis along with the clinical demands for RO water and determines the size and types of pre-treatment required to give your water system optimal performance.



Softener

Examples of Pre-Treatment Components

- Blend Valve
- City Booster Pumps
- Carbons
- Softeners
- Sediment Filtration



Pre-Treatment Wall Mount Header

Accurate pre-treatment, scheduled feed water testing, and product water monitoring, along with proper cleaning and disinfecting, are preventive measures that help to ensure the quality of the product water and extend the life of your RO system.

Solutions. Service. Support

615.355.6063

www.betterwater.com

