

BETTER WATER, INC[®]

FLOOR BOXES - MEDICAL DEVICE

FLOOR BOX



For effective containment of dialysis machine connections, Better Water, Inc has designed the Floor Valve Box. This component serves as an easily accessible container which conveniently houses all your dialysis machine connections. With the BWI Floor Valve Box your dialysis work area is kept free from leakage and corrosion.

The Better Water, Inc Floor Valve Box is made of light-weight thermoplastic. Thermoplastic is corrosion-free, durable and easily cleaned. The BWI Floor Valve Box can accommodate up to four loops for primary and secondary concentrate, bicarbonate and pure water. For easy drain connections, the unit is equipped with a drain tail piece.

Standard Features


- Easy access to machine connections
- Made of light-weight Thermoplastic
- Effectively contains leakage
- Corrosion Free
- Connections are quick connect / disconnect
- Built in check valves to minimize spillage
- Concentrate and Bicarbonate regulators are included
- Tubing, fittings and valves are included for connection to loops



Solutions. Service. Support
www.betterwater.com

Equipment changes are at the sole discretion of BWI. Brochure photographs & specs may not reflect recent changes.

FLOOR BOXES



BWI carries over 2000 parts, accessories and consumables in stock, and ready to ship including:

- Carbon Filters (provides 10 Minute Empty Bed Contact Time)
- Storage Tanks
- Duplex & Simplex Softener
- Brass pressure regulator
- Blend valve
- City boost pump
- Thermometer
- Floor valve boxes
- Remote alarms
- & More!



Central Water Systems Central Delivery Systems Reverse Osmosis Units Portable RO Units

Pre-Treatment Racks Post Treatment Racks D.I. Exchange Accessories Consumables

Hemodialysis Medical Pharmaceutical Laboratories Commercial Industrial

Better Water, Inc
 698 Swan Drive Smyrna, TN 37167
 Phone: 615.355.6063 Fax: 615.355.6065
www.betterwater.com